



Camp hip abduction
orthosis for children

Camp hip abduction orthosis for children

Intended use is to support maturing of the hip with hip dysplasia and used after re-posed hip dislocation with infants. The orthosis fixes the hip joint in a 90° flexion and 60° abduction position.

Indications

- Hip dysplasia (ICD10: Q65.8)
- Hip joint luxation (ICD10: S73.00)
- Hip joint subluxation (ICD10: T14.3)

Product function

- Limiting range of motion
- The legs are positioned to secure and help mature the hip bone and to prevent Hip (sub) Luxation

Features

- Sideways adjustable hip joints
- Foam padded Thigh Cuffs with Velcro closure
- Chest strap with Velcro closure
- Material: back part made of plastic with polyethylene foam lining
- Color: white



Sizing chart

Article number	Spread width
2901	12 - 16 cm
2902	16 - 23 cm
2903	18 - 26 cm
2904	23 - 34 cm
2905	30 - 43 cm